

# Surgery Platforms

## Surgery Platform



**USE:** Can be used for surgery, micro manipulations, general micro-utility platform (i.e. organ inspections, dissections, etc).

- Two Flexible LED lights - provides direct illumination with no shadows and no heat.
- Flex arm, non-obstructive rimless magnifying glass (main lens: 2.3 X, inset lens: 4.0 X). (Allows continuous eye contact with object when moving into magnification area).
- Limb restraints, Non-skid feet for added stability.
- Size: 9" X 13.5", Universal power supply (100 – 240 VAC, 50-60 Hz).

**ORDER:** 14370123

## Thermo Controlled Surgery Platform



- Two Flexible LED lights - provides direct illumination with no shadows and no heat.
- Flex arm, non-obstructive rimless magnifying glass (main lens: 2.3 X, inset lens: 4.0 X). (Allows continuous eye contact with object when moving into magnification area)
- 4 Space Gels for continuous heat, 2 Space Blankets to retain heat.
- Thermo Reflector & Turbo Heat Sink to maximize heat transfer.
- Limb restraints, Non-skid feet for added stability.
- Size: 13.5" (L) x 9" (W) Universal power supply (100 – 240 VAC, 50-60 Hz).

**ORDER:** 14370124

*Maintain Animal Body Temperature During Surgery !*

Ultra-Safe Thermo Regulation:

- will not over heat animal
- will not burn animal
- no electricity used
- no bulky water bath
- no water bed effect

## Thermo Recovery Platform



- Convenient size 13.5" (L) x 7.5" (w), provides over 2 hrs of heat.
- 6 Space Gels for continuous heat, 4 Space Blankets to retain heat.
- Thermo Reflector and Turbo Heat Sink to maximize heat transfer.
- Fits into a standard 18 (L) x 9 (W) shoe box cage.

**ORDER:** 14370125

*Animal Recovery Made Simple and Easy !*

## Gel Pad Reactivation

- Small size is ideal for lab bench 12" (L) x 10" (W) x 8.5 " (H).
- Easily reactivate four gel pads in a single run, (Place in chamber, start, remove pads after 15 min)
- Safe and easy to operate **ORDER:** 14370140

*Quick, Easy, Safe*



**Fisher Scientific**

For customer service, call 1-800-766-7000.  
To fax an order, use 1-800-926-1166.  
To order online: [www.fishersci.com](http://www.fishersci.com)

BN0113117