New

CO₂ Chamber

for mice and rats

Door

- Place animals in chamber without losing CO₂.
- Provides for ultra quick animal induction.



Simple, Quick & Easy to Use!

The Cleanest CO2 Induction Chamber on the Market!

Delivery Tube Within Chamber

• Fills chamber from bottom, decrease induction time, "more humane".

Quick Connect w/valve

Closes when disconnected, keeps CO₂ in tube and chamber.

7' Tube

Easily connects to CO₂ source.

Clear Acrylic Top

 Easy to see animals in chamber.



Solid Sides

 Prevents viewing of animals during the process.

Small Foot Print

 9" diameter, easy to keep on bench.





Disposable Bag Liner

- Animal feces & urine <u>remain</u> in bag with animal.
- Chamber is always clean.

Cat #: COIC CO₂ Chamber

Includes: See thru acrylic lid with door & quick connect w/ valve, 7' tubing, 50 disposable insert bags.



Cat #: COBL Disposable Bag Liner (50 bags/pk)

Use w/ CO₂ Chamber.

To ORDER: Contact your VWR sales agent! VWR.com 800-932-5000